fiocchetti

SCIENTIFIC REFRIGERATORS

Armadi frigoriferi termostatati

Thermostatic refrigerator cabinets

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
At any time, the manufacturer may bring modifications to the model above described, for technical or commercial reasons. The picture does not necessary shows the product in the standard version.

ECT-F CONTROL
DMLP DIGITAL MONITOR
DMLP PRINTER

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
*functional
safe
performing
interactive
cost-effective*

supplyLAB
www.supplylab.it

standard

ECT-F CONTROL Riduce i consumi Reduces power consumption

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Fiochetti's new refrigerators and freezers fitted with microprocessor ECT-F CONTROL guarantee excellent performances and at the same time an effective power saving, when compared with previous models.

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INTERACTIVE

Thanks to a RS485 IN/OUT port and some available accessories, our equipment may be monitored from remote through ETHERNET or GSM. Moreover, thanks to some available functional analysis from remote, it is possible to intervene on the working parameters of the control system, solving any eventual problem.

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THE NEW COOLING MANAGEMENT SYSTEM

FUNCTIONAL, thanks to the easy-accessible and friendly interface, both for end user and for Service. **SAFE**, with alarms and errors occurrences list and the special and exclusive functions developed by us. **PERFORMING**, thanks to the advanced control of defrosting phases and of cooling circuit. Our new system guarantees to our customers the maximum flexibility and reliability, keeping their products and investment safe.

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
Latest generation electronic controller, with 0,1 °C resolution, fitted with backlit alphanumeric display with a 5 languages menu.

Clear and precise information concerning the working status of the unit is displayed, and so are all the detected critical temperatures. A password, when activated, identifies two separate access level to the controller (end user and Service), thus protecting the working parameters, the On/Off main switch (on request) and the operating set-point of the unit. The interactive display may be customized and enables an immediate check of all main settings for a fast trouble shooting of eventual problems (Doctor View). Special functions such as Stop-on-door, Second cooling and Advanced defrosting functions allow an efficient managing of the defrost cycles, reducing thermal variations and consumption.

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A continuous self-monitoring on the main settings of the equipment allows to immediately identify alarms and faulty conditions. An acoustic (mutable) and a visual alarm inform the operator about the critical events due to High Temperature, Low Temperature, Door Alarm, and Power Failure (only when electric power is restored, and if the temperature has exceeded the pre-set limits).

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
The list of detected alarms, detailed with date/hour and type of event, duration and maximum temperature reached, can be easily accessed from the end user menu. Possible breakdowns of the equipment are identified and transmitted through a constant detection of the temperatures of evaporator, condenser and probes, thus reducing or completely avoiding any damage to the stored products.

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DMLP Digital Monitor is an innovative alarm system completely independent for a constant monitoring and temperature recording. DMLP Digital Monitor will work during power failure, thanks to its automatically re-charged battery, can be removed and it is fitted with a double memory (short and long term). An acoustic signal (that can be muted, and re-activated automatically) and a visual one show the user the critical events due to High temperature, Low temperature, Open Door, Power failure and fault conditions. DMLP Digital Monitor also detects the working status of the equipment, the number of open door logs and any possible modification to the working parameters. A password, when activated, identifies two separate access level to the controller (end user and Service), thus protecting the working parameters, the ON/OFF main switch, the pre-set temperature limits and the internal clock.

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Extractible **Secure Digital**
able to record the data referred
to 5 years working of the unit.

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Self-diagnostic interactive TEST

Simply pressing a button the instrument will test the system integrity by checking the state of the buffer battery, of the signalling buzzer, showing the set limits, simulating an open door alarm and verifying interactively the door micro switch.

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
DMLP Monitor is the first system ever equipped with an extractible Secure Digital card able to record the data referred to 5 years working of the unit (with 1 minute slot). SD includes pre-installed software for the download of the data on any Windows® computer. It is the easiest and less cost effective way to store working sensitive information in a digital format.

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technical data sheet

- Resolution 0,1 °C. Accuracy +/- 0,3 °C
- Back-lit alphanumeric LCD Display, with 16 characters and integrated keypad with no.4 soft-touch pads
- Menu available in 5 languages, customizable
- Turn on, turn off and set point modification accessible by password
- Buzzer
- Internal clock with independent battery
- Functions and working status with internal temperature are displayed, but it is also possible to customize the display with different visualisations mode, i.e.
 - Set point and actual temperature
 - Date and hour
 - Relative humidity and temperature (only if sensor is fitted, Optional)
- Recording of the last 12 alarms with date/hour, max reached temperature and duration of the alarm state

Alarm notification:

- **Door ajar**
- **High temperature**
- **Low temperature**
- **Power failure** (only when electric power is restored, and if the temperature has exceeded the pre-set limits)
- **Dirty condenser**

Breakdown signalling and highest level safety:

- **Faulty condenser**, in case of too high condensing temperature
- **Faulty evaporator**, in case of too low evaporating temperature

- **Faulty sensors**: in case of ambient sensor failure, the equipment will keep on working correctly until the intervention of a service engineer and with a standard working status.
- **Faulty compressor protection**: in case of malfunctioning, at first the compressor is stopped, then the controller tries to restart the unit. If the operation is unsuccessful, the equipment is definitely turned off to avoid any possible further damage.


SPECIAL FUNCTIONS (standard fitted)

- **Advanced Defrost Monitoring**: continuous monitoring of evaporator maximum and minimum temperatures during every working cycle. In case a too high ice thickness level is detected on the evaporator, the controller starts to defrost. For a better defrosting phase, the evaporator fan is shut, in order to shorten its duration and to avoid any possible temperature increase inside the equipment.
Advantages: lower power consumption and greater guarantee of temperature uniformity
- **Second cooling**: To prevent too high temperature variations inside the refrigerator, the evaporator fan is independent, and it works exploiting – when the compressor is off – the latent cold accumulated into the evaporator
Advantages: Internal stable temperature and minor compressor working, resulting in lower power consumption.
- **Stop-on-door**: This function enables to avoid an excessive decrease of the evaporator temperature due to switching off of the fan during door opening.
It affects the compressor turn off in case of a too high evaporator overcooling
Advantages: Minor thermal shock, after door opening.
- **Doctor view**: by pressing a single button the user may visualize some important information concerning the main settings of the refrigerator (set point, temperature alarm limits and open door alarm), verify the buzzer and give immediate support to the service engineer in case of malfunctioning (instant temperature of evaporator and condenser)

SERVICE INFORMATION

- Recording of the last 32 defrosting phases, 32 days with open door logs, and 32 malfunctioning events with date/hour and duration of the event.
As far as door logs, all the openings are stored, but only the critical (i.e. those that caused an Open door alarm) are identified.
- Instantaneous reset of the producer's parameters configuration. If the parameters are modified by unauthorized personnel, it is possible to return immediately to the correct configuration
- Function INFOTEST: it's a high effective test instrument to verify the functioning of the equipment.
After the test has been run, the following working information will be displayed:
 - Average time of running/switched off compressor
 - Duration of defrosting phases
 - Duration of recover phase from defrost
 - Data concerning the heat emission and cooling power
- It is also possible to access the following important data:
 - Total no. of days and hours of equipment functioning
 - Total no. of days and hours of compressor working
 - Working time percentage (%) of compressor in the last hour and in the last 24 hours
 - Average working and shut down time of compressor
 - Number of turning on from panel and from mains
 - Maximum and minimum evaporator temperature
 - Condensing temperature increase during cooling phase
 - Installed firmware release
 - Identifying serial number of compressor

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SD includes pre-installed software for the download of the data on any Windows® computer.

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DMLP Monitor:

- Resolution 0,1 °C. Accuracy +/- 0,3 °C
- Back-lit alphanumeric LCD Display, with 16 characters and integrated keypad with no.4 soft-touch pads
- Menu available in 5 languages, customizable
- Turn on, turn off and set point modification accessible by password
- Buzzer
- 128 Mb SD card
- Internal clock with independent battery
- Independent check (alarm) of the battery level and constant control of tension
- No.2 output relays (dry contact)
- RS485 IN/OUT port, for additional connectivity modules
- RS232 Port for DMLP Printer
- Port for DMLP PT100 Extension

DMLP monitor can be visualized through a remote box

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DMLP Printer:

(available only in addition to DMLP Monitor).

A modern and flexible answer to the graphic recorder with paper round charts.

It shelters from the risk of data loss (in case of missed substitution of paper chart, pen nib or battery, storage problems), and guarantees a perfect reading of the critical product situations. It performs an automatic weekly print (just as a graphic recorder), an automatic print of alarms reports after the critical condition is over (with date, hour, alarm duration and maximum temperature reached), a daily print, a continuous print to monitor specific processes, and several print on request (parameters, open door logs, alarm and faults lists).

Temperatures can be printed both in tabular and graphic modes, with different resolution degrees, customizable from user.

DMLP Printer is the most complete instrument for data storage on paper, inkless and fed by common thermographic paper for low running costs.

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PT100 Extension module:

(available only in addition to DMLP Monitor).

Maximum precision for temperature tracing.

It's an optional module to interface modern temperature control instruments (thanks to 4/20mA outputs) to a system based on PT100 4 Wires A Class.

The ideal solution for any equipment needing a continuous validation when installed with pre-set connections to the customer's unit.

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
Web Light Server and GSM Communicator:

available both with ECT-F Control and DMLP Monitor, it allows through the fitted RS485 port to monitor and to intervene on the units fitted with these accessories. Various access levels enable the user and service engineers to test and give instructions to the refrigerators using their mobile phones or through pc. A friendly user software, makes it easy to access the necessary information: your equipment, at a glance, wherever you are.

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Company with quality system certified by BUREAU VERITAS ISO 9001

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