



Fabbricatori di ghiaccio granulare - Flakers - Machines á glace en grains - Flockenesbereiter

supplyLAB

Cacém Park - Edifício 9
Estrada de Paço de Arcos nº88
2739-512 Agualva Cacém
T +(351) 21 4278700
F +(351) 21 4278709

www.supplylab.pt
geral@supplylab.pt

Caratteristiche tecniche
Technical features
Caractéristiques techniques
Technische Daten



Modello - Model - Modèle - Modell	BF80A		BF80W		AF10A		AF20A		AF20W		AF30A		AF30W			
Dimensioni (LxPxH) (mm) Dimensions (WxDxH) Dimensions (LxPxH) Abmessungen (BxTxH)	529X626X794				584x613x1006				584x613x1006				560x800x1465			
Piedini (H fabbricatore) (mm) Feet (Ice maker H) Pieds (H Producteur) Füsse (Bereiter H)	 (H=919)				 (H=1136)				 (H=1136)				 (H=1550)			
Condensazione Condensation Condensation Kühlung	Aria	Acqua	Aria	Acqua	Aria	Acqua	Aria	Acqua	Aria	Acqua	Aria	Acqua	Aria	Acqua		
	Air	Water	Air	Water	Air	Water	Air	Water	Air	Water	Air	Water	Air	Water		
	Air	Eau	Air	Eau	Air	Eau	Air	Eau	Air	Eau	Air	Eau	Air	Eau		
	Luft	Wasser	Luft	Wasser	Luft	Wasser	Luft	Wasser	Luft	Wasser	Luft	Wasser	Luft	Wasser		
Produzione ghiaccio (kg/24h) Ice production Production de glace Eisproduktion	80				85				120				200			
Capacità contenitore (kg) Bin capacity Capacité bac Lademenge	25				32				32				90			
Potenza (W) Power Puissance Leistungsaufnahme	330				400				500				760			
Alimentazione elettrica Electric voltage Alimentation électrique Netzspannung	230/1/50				230/1/50				230/1/50				230/1/50			
Gas Refrigerante Refrigerant gas Gaz réfrigérant Kühlmittel	R134A				R134A				R134A				R134A			
Acqua - Water Eau - Wasser	15°C (60°F)															

Legenda - Caption - Légende - Zeichenerklärung: □ di serie - as standard - standard - serienmäßig

● disponibile a richiesta - available upon request
disponible sur requête - auf Wunsch verfügbar

flocchetti
THE COLD MANUFACTURER

CF di Ciro Flocchetti & C. s.n.c. via Lorenzini, 53 - 42045 Luzzara (RE) Italy Tel. +39 0522 976232 - Fax +39 0522 976028 info@flocchetti.it www.flocchetti.it